

PATENTED RELIEF HANDLE
POST-FALL ADJUSTMENT SYSTEM

Achieve gravity override to Chair in the Air position in seconds, relieving the pressure on femoral arteries.


- Enhanced Comfort Backpad
- Quick Connect Chest
- Tongue Buckle Legs



LANYARD KEEPERS
Keeps the lanyard secure and out of the way when not in use



ENHANCED COMFORT BACK PAD WITH D-RING
Cushioning on the back and shoulders helps distribute pressure, improving users' comfort.



STRAP KEEPERS
Keeps excess webbing organized and out of the way



QUICK CONNECT CHEST BUCKLE
Easy connection and fast adjustment



LABEL PROTECTION
Durable cover protects and easily reveals labels



IMPACT INDICATOR
Clearly shows if the harness has been involved in a fall



LITEFIT™ PLUS FULL BODY SAFETY HARNESS

STANDARD



SMALL

H312001XQP Standard
H332001XQP Positioning
H332101XQP Construction

X-LARGE

H312004XQP Standard
H332004XQP Positioning
H332104XQP Construction

POSITIONING



MEDIUM/LARGE

H312002XQP Standard
H332002XQP Positioning
H332102XQP Construction

XX-LARGE

H312005XQP Standard
H332005XQP Positioning
H332105XQP Construction

CONSTRUCTION

